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- **Primus**
- **Manuell**
- **Primus-Stabil 3**
- **Manuell-Stabil 3**

CE 0123

Welch Allyn GmbH & Co. KG

Zollerstrasse 2-4 D-72417 Jungingen, Germany
 Telefon (++49) 7477-92 71-0 Telefax (++49) 7477-92 71-90
 E-Mail: info@welchallyn.de
 homepage: www.welchallyn.de
 www.speidel-keller.de

WelchAllyn
Speidel + Keller

Garantie

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Warranty

We provide against defect in materials and workmanship for your WA manometer movement the following guarantee:

- Primus and Manuell: 1 year guarantee/ • Primus Stabil 3 and Manuell 3: 3 years guarantee.
- Should the manometer be defective or deviate from the manufacturer's specifications, WA will repair or replace the unit or its parts free of charge. This guarantee does not apply in the case of damage caused by incorrect handling, improper use, accidents, modifications or during transportation. Guarantee claims shall not apply if the device is not used in accordance with the manufacturer's recommendations or is opened by unauthorized persons.

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- Primus
- Manuell
- Primus-Stabil 3
- Manuell-Stabil 3

Dear customer,

we are delighted that you have decided to purchase a precision blood pressure meter from Welch Allyn.

In this blood pressure meters, the manometer is directly connected to the pump valve. It is a typical modern one-hand unit. The bulb is attached to a lateral metal bracket, which stabilizes the unit and allows rapid inflation of the cuff. The PRIMUS units have a single tube system, while the MANUELL units have a double tube system. The cuffs are equipped with a Luer connecting piece for rapid connection or changing of the cuffs. The units are equipped with a precision screw valve. This valve is easy to use. By turning the knurled screw slightly, the air can be released at a finely adjustable rate. When this screw is turned farther to its final position, this guarantees rapid release of air which means that the air escapes very quickly until the cuff is completely empty. Before beginning the next measurement, the valve should again be closed so that the cuff can be pumped up.

Important:

All units are equipped with an excess pressure prevention device (Pressostabil®).

The Primus-Stabil 3 and Manuell-Stabil 3 have unique doublewall housing made of impact polymer in order to protect the measuring element of these thoroughly tested manometers from impact and bumps.

Fitting the cuff

After choosing the appropriate cuff size, the arm cuff is put in position, usually on the left upper arm. The arm circumference range is indicated on the cuff (also see table). The cuff should be tightly placed against the soft parts of the arm, but should not exert any pressure on the artery when deflated. The cuff should be placed approximately in the middle of the upper arm, in such a way that the tube points downward (in a distal direction) and is about in the middle of the bicep. The arm should be almost fully extended at the elbow and the muscles should be completely relaxed when the cuff is being put in place and during measurement.

Measurement

After putting the cuff in place, air is pumped into the entire system with the valve closed. After pumping up the cuff, the patient must remain completely motionless and one should avoid pressing on the cuff, as this affects the measurement. The release screw of the valve is then opened slightly in order to allow the air to escape slowly from the entire system, usually at a rate of 2-3 mmHg/second. This device is a so called one-hand model, leaving one hand free for auscultation of systolic and diastolic blood pressure with the stethoscope. After noting the upper and lower blood pressure values, the screw valve should be completely opened. This guarantees that when the cuff is put back in place, the air will have completely escaped. The valve should again be closed before beginning another measurement.

Specifications

Measurement range	0-300 mmHg
Precision	± 3 mmHg
Scale graduation	2 mmHg
Excess pressure prevention device	Pressostabil® (prevents pressure of more than 350 mmHg)

Cuff sizes for WA

Adults, large upper arm
Adults, upper arm
Children/adolescents
Infants/small children
Newborns
Adults thigh

Arm size for velcro cuffs

34 - 52 cm
29 - 42 cm
20 - 29 cm
14 - 20 cm
10 - 14 cm
41 - 60 cm

Range for hook cuffs

large adult
adult
child range
infant
newborn range
thigh range

The WA dimensions for arm size were developed over many years of cooperation with our customers and according to their wishes. WA cuffs offer optimum measurement reliability and user-friendliness.

Our One-Piece-Cuffs apply to the AAMI standards for suitable fitting.

Care instructions

Reusable Welch Allyn cuffs can be washed in a washing machine at 60°C, but must not be spin-dried. The blood pressure meter can be disinfected by wiping clean with conventional agents.

Important note

As a rule, blood pressure meters should be technically checked every 2 years. In such cases, return your unit to Welch Allyn or to your dealer.

CE-Certification

The conformity of this product with the requirements of the CE-guideline 93/42 EEC, council from june 14, 1993 concerning Medical Products and Medical-Product-Legislation from August 2, 1994, was established through a conformity check process. Therefore Welch Allyn introduced an entire Quality-Management-System in accordance with EN 9001 and EN 46001.

Calibration

Should the manometer fail to be accurate to ± 3 mmHg. WA will reset the manometer free of charge within a period of 1 year from the purchase date.

Accessories

we provide a 1 year guarantee against defects in materials or workmanship defects for cuffs, rubber bags, and tubes.

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